FINESSE[®]+

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Referenser

1. Wright TC, Massad LS, Dunton CJ, et al, 2006 consensus guidelines for the management of women with abnormal cervical cancer screening tests, *Am J Obstet Gynecol*, 2007; 197(4):346-55.
2. Wright TC, Massad LS, Dunton CJ, et al, 2006 consensus guidelines for the management of women with cervical intraepithelial neoplasia or adenocarcinoma in situ, *Am J Obstet Gynecol*, 2007; 197(4):340-45.
3. "2006 Consensus Guidelines for the Management of Women with Cervical Intraepithelial Neoplasia or Adenocarcinoma in situ." American Society of Colposcopy and Cervical Pathology (ASCCP). October 2007. Available at http://www.asccp.org/Portals/9/docs/pdfs/Consensus%20Guidelines/algorithms_hist_07.pdf. Retrieved June 27, 2012.
4. Electrosurgery for HPV-related diseases of the lower genital tract, Wright TC, Richart RM, Ferenczy A, Arthur Vision, Inc., New York, 1992.
5. Large loop excision of the transformation zone, Prendiville W, Chapman & Hall, 1993.
6. Baggish MS, Ferenczy A, Gerbie M, Richart RM, Potential complications of loop excision, *Contemp OB/GYN* 1994; 39(8):93-107.

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Technical Data

Physical Specifications

Dimensions: 14.0" (35.6 cm) W x 14.7" (37.3 cm) D x 7.3" (18.5 cm) H
 Weight: 24 lbs. (11 kg)

Regulatory Information

IEC classification: Class I, Type BF 
 Defibrillator protected

Patient leads: RF Isolated 

Mode of Operation: Continuous operation with intermittent loading (10s/30s)

Ingress protection: IPxx (ingress protection rating is not required)

Protected against spillage per IEC 60601-1 and IEC 60601-2-2

EMC: Complies with IEC 60601-1-2:2007 and IEC 60601-2-2:2009 requirements for electromagnetic compatibility (EMC) of medical electrical equipment. For detailed information on EMC, see the Electromagnetic Compatibility section of this manual.

Supply Voltage and Current Considerations

The FINESSE+ system can be obtained in 115 VAC (item no. FIN-110) or 230 VAC (item no. FIN-220) configurations.

Supply Voltages

	<u>FIN-110</u>	<u>FIN-220</u>
Nominal Voltage:	115 volts rms	230 volts rms
Regulated Range:	103-127 volts rms	206-253 volts rms
Frequency:	50/60 Hz	50/60 Hz

Output Power vs. Supply Voltage

Within the Regulated Range listed above, output power into a 500Ω load remains ±15% or ±5W, whichever is greater; of the power delivered at the center of this range in cut and blend modes and ±30% or ±10W, whichever is greater; in coagulation mode.

Maximum Supply Current and Power

<u>FIN-110</u>	<u>System Total</u>	<u>ES Module</u>	<u>Smoke Evacuator</u>
Idle State:	0.4 amp, 50 watts	0.4 amp, 50 watts	negligible
Cut Mode:	4.1 amps, 470 watts	2.1 amps, 240 watts	2.0 amps, 230 watts
Coag Mode:	3.5 amps, 405 watts	1.5 amps, 175 watts	2.0 amps, 230 watts

FIN-220	System Total	ES Module	Smoke Evacuator
Idle State:	0.2 amp, 50 watts	0.2 amp, 50 watts	negligible
Cut Mode:	2.6 amps, 585 watts	1.1 amps, 240 watts	1.5 amps, 345 watts
Coag Mode:	2.3 amps, 520 watts	0.8 amps, 175 watts	1.5 amps, 345 watts

Overcurrent Protection

The FINESSE+ is protected by two sets of 5mm x 20mm, time-lag, low breaking capacity, 250VAC fuses (all fuses are type TxAL250V, where x is the rating capacity in Amps):
Electrosurgical System: a pair of fuses protect the leads connected to the primary winding of this module's main power transformer.

FIN-110: 2.5 amp FIN-220 1.25 amp

Smoke Evacuation System: a pair of fuses protect the smoke evacuation motor:

FIN-110: 3.15 amp FIN-220 2.5 amp

Output Characteristics

Output Frequencies

Output Frequency (all modes): 450 kHz ± 50 kHz.
Waveform Pattern Repetition Rate: 28.1 kHz ± 3 kHz.

Output Values at Maximum Settings

<i>Operating Mode</i>	<i>Output Setting^a</i>	<i>Duty Cycle^b</i>	<i>MOV^d</i>	<i>Crest Factor^e</i>
Cut	99	100%	1080	1.8
Blend 1	99	62.5%	1280	2.2
Blend 2	99	50%	1420	2.5
Blend 3	99	37.5%	1500	2.8
Coag	75	c	2400	5.8

a Output Setting specified is power in watts, plus or minus 15%, delivered into a 500Ω patient load.
b Duty Cycle is ratio of burst duration to burst-plus-rest duration.
c Coag pulse consists of two high voltage cycles followed by lower amplitude ringout for about 10 μsec, repeated at 28.1 kHz.
d Maximum Output Voltage; voltage is zero-to-peak, open circuit. Lower values are permitted. Use only accessories with a Rated Accessory Voltage (RAV) that is greater than the MOV.
e Crest factor is defined as the ratio of peak voltage to RMS ("average") voltage, and is a general quantification of an electrosurgical waveform's degree of coagulation. Data shown are typical values at 500Ω.

Table 1. Output Characteristics

Output Power vs. Load

Output display is calibrated to be relative power in watts at a load of 500Ω, ±15% or ±5W, whichever is greater.

In the pure cut and blend modes, the negative-feedback output stage attempts to maintain the output voltage under low resistance loads so that the cutting characteristics are quite uniform throughout the entire cut. Cut mode maximum output is 140 watts at 250Ω patient load. See Figure 6 for output power vs. load data.

In Coag mode the negative-feedback stage is not used in order to achieve high open-circuit voltage, which is desirable. Coag mode maximum output is 100 watts at 250Ω patient load. See Figure 7 for output power vs. load data.


Figure 6. Typical output power vs. load resistance for cut and blend modes


Figure 7. Typical output power vs. load resistance for coag mode


Figure 8. Maximum Output Voltage vs. displayed settings

Maximum Output Voltage vs. Output Setting and Mode

Figure 8 shows the Maximum Output Voltage (MOV) for each mode and output setting. Only use active electrodes and pencils that have a Rated Accessory Voltage (RAV) greater than the MOV for the selected output setting and mode.

Output Power vs. Displayed Setting

The output of the FINESSE+ system increases linearly with adjustment of the front panel output settings. Figure 9 shows the typical output power of each of the five modes over the full control range. All curves shown indicate output at 500Ω load resistance.


Figure 9. Typical output power vs. display setting for cut and blend modes and coag mode

Output Error Safety Circuit Specifications

The FINESSE+ system is equipped with a safety circuit that monitors the output signal levels and disables all system functions when the output exceeds the expected "nominal" output. The error indication can only be reset by turning off the main power switch.

The output error indicator will activate if the output power exceeds the displayed output, as defined in standard IEC 60601-2-2, Particular Requirements for the Basic Safety and Essential Performance of High Frequency Surgical Equipment and High Frequency Surgical Accessories. This threshold is strictly defined, but in the FINESSE+ it is generally characterized by a threshold that is 25W-30W above the selected output setting.

CQM Circuit Specifications

CAUTION: CQM circuits used in electrosurgical generators are an excellent tool to minimize the risk of a patient injury during electrosurgery. CQM is effective at detecting changes in pad contact to the patient, and can signal the user to a potentially unsafe condition. However, these circuits should not be relied on as a method for detecting improper pad application or as a substitute for user attentiveness during the procedure.

Always apply the dispersive pad as indicated in the dispersive pad instructions for use, and periodically make a visual check the pad application to the patient.

The FINESSE+ system is capable of monitoring dispersive pad contact to the patient. To enable this function, a split-style dispersive pad is required (see compatibility list). If a split pad is not detected, the use of a solid pad is assumed and the FINESSE+ system will not detect compromise in pad contact to the patient.

For best operation, the dispersive pad should be attached to the patient prior to plugging the pad into the dispersive pad receptacle. Once the CQM system detects the pad's contact quality to the patient, a threshold is internally set that will activate the CQM error indicator if the dispersive pad partially peels away from the patient.

When the CQM-compatible dispersive pad is properly and fully applied, tests of the FINESSE+ CQM system found that the circuit will display the error indicator when approximately 30% of the dispersive pad's surface area has separated from the patient. At the threshold of the error condition, tests performed per standard IEC/EN 60601-2-2 sec 201.15.101.5 found the following pads caused less than 6°C temperature rise on the patient's skin.

CQM-compatible pads (UTMD item numbers and descriptions):

- ES-1179 Dispersive Pad, CQM (split), precorded
- ES-1180 Dispersive Pad, CQM (split), non-corded (requires reusable cable)

WARNING: To ensure safe contact of the dispersive pad to the patient, use only pads listed as compatible with the FINESSE+'S CQM system. The use of contact monitoring pads that are not on this list may not properly signal an error condition and result in patient injury.

CQM circuit frequency: 39.5 kHz

CQM pad contact range: 10-130Ω, ±5Ω

Solid pad contact range: <10Ω±5Ω, but not greater than minimum of CQM contact range

Audio Tone Specifications

All audio tone specifications are measured in decibels, A-weighted scale, at a distance of 1m from the front of the FINESSE+ system, and are nominal values.

- Cut Mode Tone: 45-70 dBA, adjustable, 625 Hz continuous tone while activated
- Coag Mode Tone: 45-70 dBA, adjustable, 1.25 kHz continuous tone while activated
- Error Indicator Tones: 70 dBA, not adjustable:
 - CQM error: intermittent 2.5 kHz tone, 0.25 sec on/0.25 sec off
 - Pad disconnect error: intermittent 2.5 kHz tone, 0.25 sec on/0.25 sec off
 - Cross-key error: continuous 2.5 kHz tone
 - Output safety error: continuous 2.5 kHz tone

Smoke Evacuator System

Air Flow vs. Supply Voltage

At the Nominal Supply Voltage, the air flow through a new disposable filter coupled to a 30 cm long 6 mm ID tube is not less than 100 liters per minute (3.5 cubic feet per minute) at the "high" setting and not less than 70 liters per minute (2.5 cfm) at the "normal" setting.

Smoke Evacuation Duration

The smoke evacuation system will begin running immediately on activation of the electrosurgery module and remain running five seconds after the output power is deactivated.

Disposable Filter Cartridge

The FINESSE Filter Pack, catalog number ESU-501, consists of a pleated HEPA paper filter followed by a compartment containing activated charcoal.

The FINESSE Filter Pack, to ensure total elimination of odors between procedures, should be replaced on a daily basis. However, the three-stage filtration system has been tested to effectively remove odors and particles for up to 15 electrosurgical procedures. Over a period of time the external disposable filter pack can be a source of odor and possible viral contamination. Therefore, it is recommended that the external filter pack be changed every day or after 15 procedures if more than 15 procedures are performed in a single day.

Third-Stage Internal Filter

The FINESSE+ Internal Filter, item number SSE-500, is a pleated ULPA filter element. Annual replacement is recommended.

Particle Removal Efficiency

Spherical particles with a mean diameter of 0.1 microns were removed with a minimum efficiency of 99.999%, as determined by laboratory tests via a latex particle challenge.

Vacuum Motor Exhaust Filter

The FINESSE+ smoke evacuation motor is equipped with a HEPA grade particulate filter. This filter protects the FINESSE+'s internal electronics from any particulate that may be generated by the motor. This filter should never need cleaning or changing. However, if the smoke evacuation system flow seems degraded after replacing all other filtration components, contact Utah Medical Products for information.

Environmental Specifications

Operational Environment

Temperature:	50°F to 104°F (10°C to 40°C)
Humidity:	10% to 93%
Pressure:	700 hPa to 1060 hPa (10.15 PSI to 15.36 PSI)

Storage and Transport Environment

Temperature:	-40°F to 158°F (-40°C to 70°C)
Humidity:	10% to 93%
Pressure:	500 hPa to 1060 hPa (7.25 PSI to 15.36 PSI)

After transport and/or storage at conditions outside of the operating environment range, allow time for the FINESSE+ system to reach the operating environment — typically one hour is sufficient.

Accessories

The following accessories are approved by Utah Medical Products for use with the FINESSE+ Electrosurgical and Smoke Evacuation System. All items are available from Utah Medical Products or its authorized distributor. Contact Utah Medical Products for a complete and current listing of supplies for the FINESSE+ system.

Electrodes

Item No. Description

Conization Electrodes^{ab}

- DCE-110 C-LETZ Conization Electrode, Large, 23mm D x 15mm Radius (blue)
- DCE-115 C-LETZ Conization Electrode, Medium, 18mm D x 11mm Radius (green)
- DCE-120 C-LETZ Conization Electrode, Ectocervical, 10mm D x 12mm Radius (tan)
- DCE-125 C-LETZ Conization Electrode, Small, 13mm D x 9mm Radius (yellow)

LETZ® Electrodes (11 cm shaft)^{ab}

- DLP-B11 UtahLoop® Electrode, 25mm x 8mm, with Safe-T-Gauge® (tan)
- DLP-E11 UtahLoop® Electrode, 25mm x 22mm, with Safe-T-Gauge® (canary)
- DLP-L11 UtahLoop® Electrode, 20mm x 15mm, with Safe-T-Gauge® (blue)
- DLP-W11 UtahLoop® Electrode, 20mm x 12mm, with Safe-T-Gauge® (white)
- DLP-M11 UtahLoop® Electrode, 15mm x 12mm, with Safe-T-Gauge® (green)
- DLP-S11 UtahLoop® Electrode, 10mm x 10mm, with Safe-T-Gauge® (yellow)
- DLP-SQ1 UtahLoop® Electrode, 10mm x 10mm square, with Safe-T-Gauge® (orange)
- DLP-SQ2 UtahLoop® Electrode, 10mm x 4mm, 15° bend, 11 cm shaft (brown)
- DLP-T11 UtahLoop® Electrode, 5mm x 5mm, (purple)
- DBL-511 UtahBall® Electrode, 5mm diameter (red)
- DBL-311 UtahBall® Electrode, 3mm diameter (black)
- DLP-N11 Macro Needle Electrode, 0.8 mm dia x 15mm long, 11 cm shaft (red)
- DLP-U11 Micro Needle Electrode, 0.2 mm dia x 15mm long, 11 cm shaft (black)

Extender Electrode^{ab}

- DXT-S06 DXTender® Electrode Extender, small displacement (blue)
- DXT-L09 DXTender® Electrode Extender, large displacement (blue)
- DLP-X10 Electrode Extender, 10cm (blue)

External Lesion Electrodes (5.5cm shaft)^{ab}

- DLP-B05 Round Loop Electrode, 25mm x 8mm (tan)
- DLP-M05 Round Loop Electrode, 15mm x 12mm (green)
- DLP-S05 Round Loop Electrode, 10mm x 10mm (yellow)
- DLP-T05 Round Loop Electrode, 5mm x 5mm (purple)
- DLP-SQ3 Square Loop Electrode, 10mm x 4mm, 15° bend (pink)

<i>Item No.</i>	<i>Description</i>
DLP-N05	Macro Needle Electrode, 0.8 mm dia × 5mm long (red)
DLP-U05	Micro Needle Electrode, 0.2 mm dia × 5mm long (black)
DBL-505	Ball Electrode, 5mm diameter (red)
Blade Electrodes^{ab}	
CBE-100	Epitome® .4 Scalpel Electrode, standard length (blue)
CBE-150	Epitome® .4 Scalpel Electrode, extended length (blue)
CBE-200	Epitome® .2 Scalpel Electrode, standard length (white)
CBE-250	Epitome® .2 Scalpel Electrode, extended length (white)
CBE-220	Epitome® .2 Scalpel Electrode, standard length, with ZapGuard™ (black)
CBE-260	Epitome® .2 Scalpel Electrode, 4" extended length, with ZapGuard™ (black)
CBE-270	Epitome® .2 Scalpel Electrode, 6" extended length, with ZapGuard™ (black)

Microneedle Electrodes^{ab}

DN-0200	OptiMicro™ Needle Electrode, 2cm shaft
DN-0300	OptiMicro™ Needle Electrode, 3cm shaft
DN-0400	OptiMicro™ Needle Electrode, 4cm shaft
DN-0245	OptiMicro™ Needle Electrode, 2cm shaft with 3mm long 45° tip
DN-0345	OptiMicro™ Needle Electrode, 3cm shaft with 3mm long 45° tip
DN-0445	OptiMicro™ Needle Electrode, 3cm shaft with 10mm long 45° tip
DN-0800	OptiMicro™ Needle Electrode, 8cm shaft
DN-0810	OptiMicro™ Needle Electrode, 8cm shaft with 10mm exposed tip

Accessories

<i>Item No.</i>	<i>Description</i>
ESU-170	Dual Pedal Waterproof Footswitch
ESU-305 ^{ab}	Two Button Hand-switching Pen
ES-9135-LP ^b	Dispersive electrode pad with Safety Ring, solid
ES-1179 ^b	Dispersive Pad, CQM (split), precorded
ES-1180 ^b	Dispersive Pad, CQM (split), non-corded (requires ES-21174 reusable cable)
ES-21174	Reusable Dispersive Pad Cord for ES-1180 non-corded dispersive pad
ESU-501 ^c	FINESSE Filter Pack (includes speculum tubing)
ESU-502 ^b	Flexible Speculum Tubing with adapter
951-712	Universal Smoke Evacuation Tubing Set
SSE-500	FINESSE+ Internal Filter
SSE-610	Internal Filter Retaining Ring
0007-503	Fuse, 3.15 amp, 5mm × 20mm, Time Lag, for FIN-110 smoke evacuation module (type T3.15AL250V)
0007-502	Fuse, 2.5 amp 5mm × 20mm, Time Lag, for FIN-110 electrosurgical module and FIN-220 smoke evacuation module (type T2.5AL250V)
0007-501	Fuse, 1.25 amp 5mm × 20mm, Time Lag, for FIN-220 electrosurgical module (type T1.25AL250V)

- a supplied sterile
- b disposable, intended for single patient use only
- c disposable, rated for 15 procedures

Glossary

This glossary defines words and expressions used in this manual which may be new to the beginning practitioner of electrosurgery, as well as more common words that have an unusual meaning in the art of electrosurgery.

Activated charcoal: Charcoal granules that have been heated in steam for a sustained period of time so that the surfaces are free of organic gas molecules. These granules have a very large capacity for adsorbing organic gases that contact their surfaces.

Active electrode: The electrode used as a tool at the surgical site, e.g. the loop.

Adsorption: The assimilation of gas, vapor, or dissolved matter by the surface of a solid or liquid. See also activated charcoal.

Amp: Short for Ampere. The measure of electron flow in an electrical circuit. One amp represents 6.24×10^{18} electrons passing a given point in one second.

Blend mode: A mode of electrosurgical cutting wherein a periodically interrupted voltage waveform is applied to the active electrode, which causes sufficient heating of the surface exposed by the cut that desiccation coagulation occurs.

Capacitive coupling: Effective connection from one circuit to another by capacitance rather than conductance. Capacitance between circuits may accidentally come about when a conductive surface connected to one circuit lies intimately close to a conductive surface connected to the other circuit, even though the two conductive surfaces are separated by an insulating material or air. In this situation, high frequency current may flow from one circuit to the other.

Capacitive grounding: Capacitive coupling to ground. For example, in electrosurgery, the patient dispersive electrode, which was designed to be isolated from ground, may be effectively grounded by falling off the patient and lying flat on a metal table. Even though some intervening insulating material may prevent the conductive surface of the pad from actually touching the table there is nevertheless capacitive coupling and current can flow between the pad and the table. If another part of the table should be touching the patient while cutting or coagulation is being performed, a patient burn could occur at the point of contact.

Carrier frequency: The frequency of the principal radio frequency waveform produced by the electrosurgical generator.

Coag/Coagulation: The process of causing bleeding to cease by thermally clotting blood and shrinking vessel ends. Also the setting on an electrosurgical generator that produces a waveform that is effective in bringing about coagulation.

Contact Quality Monitoring (CQM): A safety circuit used in many electrosurgical systems designed to detect partial detachment of the dispersive pad from the patient's skin. When using a non-CQM dispersive pad (typically these have a single conductive surface), pad detachment can

contribute to a potentially unsafe condition, since the electrosurgical return current will be more concentrated (i.e., have a higher current density) and result in an increase of skin warming. Continuing to operate an electrosurgical generator with a partially detached dispersive pad has caused severe burns at the dispersive pad site. CQM circuits have shown to be an effective means of reducing the incidence of pad site burns due to pad detachment. Standard IEC 60601-2-2 Particular Requirements for the Safety of High Frequency Surgical Systems addresses the requirements for CQM system performance.

Current density: The degree of concentration of electrical current through a conductive surface or cross section. High current densities cause high power densities, causing high heat densities, causing high temperatures.

Cut: Separation of tissue. In electrosurgery, a process wherein a narrow line of cells are heated to the point of bursting as a thin active electrode or loop is moved through the tissue. Also the setting on the electrosurgical generator that will produce the best waveform for cutting with minimal coincident coagulation.

Desiccation: Drying. In electrosurgery, a mode of coagulation wherein cells are electrically heated in the vicinity of an active electrode to the point that they are dried and shrunk.

Dispersive electrode: Also referred to as dispersive pad. An electrode applied to the patient's skin, having a large area of application so that the current flowing between the patient and the electrode is diffuse rather than dense. A properly applied and connected dispersive pad is the most important consideration in preventing burns in electrosurgery.

Duty cycle: A parameter in the description of a repetitive waveform which gives the percentage of time that voltage is applied to the circuit compared with the total time over which a cycle occurs.

Electromagnetic Compatibility: The ability of electrical equipment to operate satisfactorily in the presence of other electrical equipment ("electromagnetic immunity"), without introducing radiofrequency interference that would affect that other equipment ("electromagnetic emissions").

Frequency: As pertaining to an electrical waveform, the rate of repetition of identical voltage or current patterns. Frequency is measured in Hertz (Hz), cycles per second.

Fulguration: Flashing, sparking. In electrosurgery, a process wherein sparks are caused to jump across a gap between the active electrode and the patient tissue where coagulation is desired.

Isolated leads: Electrical leads upon which the voltages carried are not referenced to earth ground. Perfectly isolated leads, even though they may carry large currents at high voltages, will not convey current to earth grounded conductors that may come in contact with them or the patient to which they are connected.

Leakage current: Current flowing in a path other than the intended path. For this discussion, currents that may flow through the patient or through the bodies of attending care personnel. Depending upon their magnitudes and frequencies, leakage currents may cause shocks, burns, or fibrillation.

Load resistance: The resistance to electric current between the active electrode and the dispersive or return electrode. In actual electrosurgery, the resistance presented to current passing through the patient. In test conditions it is the resistance of a resistor of known value placed across the patient leads from the generator.

Maximum Output Voltage: The maximum possible peak voltage between patient circuit connections. For monopolar electrosurgery systems, this is typically the peak voltage between the active electrode and the dispersive electrode.

Modulation: Variation in the amplitude of a waveform. Modulation of a continuous cutting waveform creates a blend waveform that will provide shallow coagulation as well as cutting. Greater than 100% modulation creates a blend mode waveform consisting of a sequence of repeated pulses each containing a number of cycles of the carrier frequency.

Monopolar electrosurgery: Electrosurgery wherein only a single active electrode is used at the surgical site, as opposed to bipolar electrosurgery wherein current flows between two narrowly separated active electrodes at the surgical site.

Ohm: The measure of electrical resistance. One ohm is the amount of resistance that requires one volt to cause a current of one amp. Represented as Ω .

Radio frequency: Frequencies of alternating electromagnetic currents higher than 20 kHz.

Rated Accessory Voltage: The maximum peak voltage that can be applied to the active electrode, relative to the dispersive electrode. The use of an accessory where the Rated Accessory Voltage is less than the Maximum Output Voltage may lead to an unsafe operating situation and patient and/or user injury.

Sinusoid Waveform: A smoothly-oscillating waveform. A continuous sinusoid is characteristic of an electrosurgical "pure cut" mode, and an interrupted sinusoid waveform is characteristic of a blended cut waveform.

Spark: When a very high voltage is applied across a short gap containing a gas (e.g. air or water vapor), the gas molecules in the gap are ionized and current flows. When this happens, light, sound, and heat are produced. This event is a spark, also called an arc.

Spray coagulation or spray fulguration: A method of electrosurgical coagulation wherein a newly exposed surface is treated by spraying it with sparks from an active electrode. In this method the electrode is moved parallel to the tissue surface, but separated by a slight distance so that the electrode is not touching but is close enough for sparks to jump. The waveform used has a very high voltage and a very short duty cycle.

Volt: The measure of electrical potential or force. Higher voltage causes higher current to flow through a given resistance. Very high voltage can cause current to flow through a gas or air gap by ionizing the gas molecules in the path.

Watt: The measure of electrical power. Power in watts equates directly to rate of heat production. Concentration of power means concentrated heating which brings about high temperature.

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
Electromagnetic Compatibility

The FINESSE+ system has been designed to operate with minimal interference to nearby electronics (electromagnetic emissions), and also to perform without any effects in the presence of other electronics (electromagnetic immunity). Because electrosurgical systems use high frequency energy for therapeutic surgical effects, the following tables are provided to satisfy regulatory requirements, and so that users of the FINESSE+ have guidance on operating the system in an environment containing multiple electronic systems.

Table 3. Guidance and manufacturer's declaration — electromagnetic emissions		
The FINESSE+ System is intended for use in the electromagnetic environment specified below. The customer or the user of the FINESSE+ System should assure that it is used in such an environment.		
Emissions Test	Compliance	Electromagnetic environment — guidance
RF emissions CISPR 11	Group 2	The FINESSE+ System must emit electromagnetic energy in order to perform its intended function. Nearby electronic equipment may be affected.
RF emissions CISPR 11	Class A	The FINESSE+ System is suitable for use in all establishments other than domestic, and may be used in domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes, provided that the following warning is heeded: Warning: This equipment/system is intended for use by healthcare professionals only. This equipment/system may cause radio interference or may disrupt the operation of nearby equipment. It may be necessary to take mitigation measures, such as re-orienting or relocating the FINESSE+ or shielding the location.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Table 4. Guidance and manufacturer's declaration — electromagnetic immunity			
The FINESSE+ System is intended for use in the electromagnetic environment specified below. The customer or the user of the FINESSE+ System should assure that it is used in such an environment.			
Immunity Test	IEC 60601 test level	Compliance level	Electromagnetic environment — guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV line(s) to line(s) ±2 kV line(s) to earth	±1 kV line(s) to line(s) ±2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.

Table 4. Guidance and manufacturer's declaration — electromagnetic immunity (cont.)

Power frequency (50/60Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in U_T) for 0.5 cycles 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 s	<5% U_T (>95% dip in U_T) for 0.5 cycles 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the Finesse+ System requires continued operation during power mains interruptions, it is recommended that the Finesse+ System be powered from an uninterruptible power supply or battery.
NOTE: U_T is the a.c. mains voltage prior to application of the test level.			
Conducted RF IEC 61000-4-6	$3V_{rms}$ 150 kHz to 80 MHz	$3V_{rms}$	Portable and mobile RF communications equipment should be used no closer to any part of the FINESSE+ System, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter: Recommended separation distance $d = 1.2\sqrt{P}$ $d = 1.2\sqrt{P}$ 80 MHz to 800 MHz $d = 2.3\sqrt{P}$ 800 MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey, ^a should be less than the compliance level in each frequency range. ^b Interference may occur in the vicinity of equipment marked with the following symbol: 
Radiated RF IEC 61000-4-3	3 V/m 80 MHz to 2.5 GHz	3 V/m	
NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the FINESSE+ System is used exceeds the applicable RF compliance level above, the FINESSE+ System should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the FINESSE+ System. b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3V/m.			

Portable and mobile communications equipment can affect electromedical equipment. Table 5 provides recommended separation distances for the operation of these devices in the vicinity of the FINESSE+ system. These immunity distances were established from the testing procedures and guidance of international standard IEC 61000-4-3 and IEC 61000-4-6.

<i>Table 5. Recommended separation distances between portable and mobile RF communications equipment and the FINESSE+ System</i>			
The FINESSE+ System is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the FINESSE+ System can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the FINESSE+ System as recommended below, according to the maximum output power of the communications equipment.			
Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter		
	150 kHz to 80 MHz $d = 1.2\sqrt{P}$	80 MHz to 800 MHz $d = 1.2\sqrt{P}$	800 MHz to 2.5 GHz $d = 2.3\sqrt{P}$
0.01	0.12 m	0.12 m	0.23 m
0.1	0.38 m	0.38 m	0.73 m
1	1.2 m	1.2 m	2.3 m
10	3.8 m	3.8 m	7.3 m
100	12 m	12 m	23 m
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter; where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer: NOTE 1: At 80 MHz and 800 MHz, the higher frequency range applies. NOTE 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			

Examples of approximate frequencies and power outputs for common RF communications devices:

- Mobile Phone: 0.01 to 1 W 850 MHz to 2.1 GHz
- Wireless Router: 0.1 W 2.4 GHz
- Cordless phone: 0.01 W 1.9 GHz to 5.6 GHz
- Laptop WiFi: 0.01 to 0.05 W 2.4 GHz
- Bluetooth Device: 0.001 to 0.1 W 2.4 GHz to 2.5 GHz

Bibliography

The following references cite some of the published literature pertaining to electrosurgical procedures. These papers should provide an introduction to the technique. Further instructions should be obtained from qualified experts before performing these procedures.

Baggish MS, Campion MS, Ferenczy A, Richart RM, *Exploring excision technique with LEEP. (panel discussion)* Contemp Ob/Gyn 37(sp):81-100 (1992)

Baggish MS, Ferenczy A, Gerbie MV, Richart RM, *Potential complications of loop excision (panel discussion)* Contemp Ob/Gyn 39 (8):93-107 (1994)

Cartier R, Sopena B, Cartier I, *Use of the diathermy loop in the diagnosis and treatment of lesions of the uterine cervix (abstract)*. Proceedings of the International Federation for Cervical Pathology and Colposcopy, London (1981)

Hainer BL, *Fundamentals of electrosurgery*, J Am Board Fam Prac; 4:419-26 (1991)

Luesley DM, Cullimore PC, Redman CW, et al, *Loop diathermy excision of the cervical transformation zone in patients with abnormal cervical smears*. Br Med J 300:1690-3 (1990)

McLucas B, Emens M, Hamou J, Rothburg R, *Diathermy loop treatment of CIN, international perspectives. (panel discussion)* The Female Pat. 15:79-90 (1990)

Mor-Yosef S, Lopes A, Pearson S, Monaghan JM, *Loop diathermy cone biopsy*. Obstet Gynecol 75:884-886 (1990)

Prendiville W, Cullimore J, Norman S, *Large loop excision of the transformation zone (LLETZ). A new method of management for women with cervical intraepithelial neoplasia*. Br J Obstet Gynaecol 96:1054-1060 (1989)

Large loop excision of the transformation zone, a practical guide to LLETZ, ed. Prendiville, W, 1st edition 1993, Chapman & Hall Medical, London.

Richart RM, *Ways of using LEEP for external lesions*, Contemp Ob/Gyn: 5:138-52 (1992)

Office Dermatology, Roenigk, Henry H., Jr., 1981, Williams & Wilkins, Baltimore

Office Practice of Skin Surgery, Schultz, Bryan C. and McKinney, Peter, 1985, W.B. Saunders Co., Philadelphia

Scoular A, *Choosing equipment for treating genital warts in genitourinary medicine clinics*, Genotourin Med; 67:413-19 (1991)

Cutaneous Electrosurgery, Sebben, Jack E., 1989, Year Book Publishers, Inc., Chicago.

Plastic Surgery, Smith, James W. and Aston, Sherrell J., fourth edition 1991, Little, Brown and Company, Boston

Soderstrom R, *Taking a colposcopic-directed loop biopsy*. Contemp Ob/Gyn 37(sp):77-78 (1992)

Wright TC, Massad LS, Dunton CJ, et al, *2006 consensus guidelines for the management of women with abnormal cervical cancer screening tests*, Am J Obstet Gynecol, 197(4):346-55 (2007)

Wright TC, Massad LS, Dunton CJ, et al, *2006 consensus guidelines for the management of women with cervical intraepithelial neoplasia or adenocarcinoma in situ*, Am J Obstet Gynecol, 197(4):340-45 (2007)

Electrosurgery for HPV-related diseases of the lower genital tract, Wright TC, Richart RM, Ferenczy A, 1992, Arthur Vision, New York.

Wyre HW, Stolar R, *Extirpation of warts by a loop electrode and cutting current*, J Dermatol Surg Oncol 3(5):520-22 (1977)

The following standard is globally recognized for safety and performance certification of electrosurgical systems, and can be referenced for the basis of many specifications. This standard is not the exclusive requirement for electrosurgical system safety, as many other standards exist for electromedical equipment and supplies. Nor should it be assumed that an electrosurgical system meets these requirements, unless the proper regulatory tests have been performed and regulatory markings are provided on the system.

Medical electrical equipment — part 2-2: Particular requirements for the basic safety and essential performance of high frequency surgical equipment and high frequency surgical accessories, IEC 60601-2-2 Edition 5.0 (2009), International Electrotechnical Commission, Geneva, Switzerland.

Warranty and Service

FINESSE+ Limited Warranty

Utah Medical Products, Inc. (UTMD) warrants each new FINESSE+ (the product) against defects in materials and/or workmanship for a period of two years from the date of purchase and agrees to repair or replace any defective product without charge. This warranty does not cover damage resulting from accident, misuse, lack of reasonable care or proper maintenance. This warranty shall be void if the product is repaired by anyone other than UTMD or an authorized service agent. This warranty does not extend to anyone other than the original purchaser nor to accessories manufactured by other vendors.

Except as provided herein, UTMD makes no warranties of any kind, either expressed or implied and specifically excluding any warranty of merchantability or warranty of fitness for a particular purpose.

UTMD will not be liable for any special, nonconsequential or incidental damages arising out of the use or inability to use the product. In no event shall UTMD's liability hereunder exceed the purchase price of the product. This warranty shall be void and of no force and effect with respect to any product which is damaged as a result of A) neglect, alteration, electric current fluctuation or accident; B) improper use, including failure to follow proper operating and maintenance instructions, and to provide proper environmental conditions prescribed in UTMD's product instruction manual; C) repair by other than UTMD or authorized service agents appointed by UTMD's service announcements; or D) use of supplies or parts which do not meet UTMD specifications.

Obtaining Service

To obtain service, please call UTMD at 800-533-4984 (outside USA and Canada – 801-566-1200) to receive specific instructions. Please be prepared to provide:

- Your UTMD account number, if known
- Your billing address
- Your shipping address
- A description of the problem with, or the service requested for, the FINESSE+ system
- A purchase order number or authorization to bill the service, if the FINESSE+ system is not under warranty.

A return goods authorization (RGA) number will be issued, and should appear on the return address label. Returns without an RGA number may be refused. To avoid shipping damage, the FINESSE+ system should only be shipped in its original packaging. If you have not retained the original packaging, UTMD can arrange for a replacement box and packing material.

Index

A

Accessories
list of 35
selection of 3
warnings 3
Activated charcoal 37
Activation of smoke evacuation system 17
Active electrode 15, 37
warnings 3
Adsorption 37
Alarms. See Errors
Alternate current paths 16
Amp 37
Annual maintenance 26
Audible tones
volume control 21
Audio tones
specifications 33

B

Ball electrode 35
Bibliography 47
Blade electrodes 36
Blend mode 10, 37
selector switch 18
Burns
electrosurgical 16

C

Capacitive coupling 37
Capacitive grounding 37
Cardiac pacemakers 4
Carrier frequency 37
Charcoal
activated 37
Cleaning 26
C-LETZ electrodes 35
Coag/Coagulation 37
Coag mode controls 19

Coag output control knob 19
Coagulation 15
desiccation 10, 38
fulguration 39
fulguration vs desiccation 10
spray 39
Complications
electrosurgical 16
Connectors
dispersive pad 20
front panel 20
handswitch 20
rear panel 21
Connectors and indicators
front panel 17
Contact Quality Monitoring (CQM) 38
Continuous sinusoid 39
Controlled Output Circuitry 10
CQM
error 8, 20
procedure for using 24
Crest factor
specifications 30
Cross-key 7
error 19
Current density 38
Cut 38
Cut/Blend mode controls 18
Cut/blend mode select switch 18
Cut/blend output control knob 18
Cut mode
vs. blend modes 10
Cutting 15
troubleshooting 27

D

Daily maintenance 26
Desiccation 10, 38
Dispersive electrode. See Dispersive pad
Dispersive pad 38
circuit specifications 32
error 8, 19
pad type indicators 20
placement 12, 24
receptacle 20
solid 36
troubleshooting 27
warnings 2
Duty cycle 38
specifications 30

E

Electric shock 16
Electrode
active 15, 37
selection 12
Electrodes 35
Electromagnetic Compatibility 4, 38
using compatible accessories 3
Electromagnetic Interference 4, 43
from RF communication devices 45
Electrosurgery
complications 16
monopolar 15
principles 15
waveforms 16
Electrosurgical generator
design 7

- Environmental specifications 34
- Epitome Scalpel electrodes 36
- Errors 7
 - CQM 8, 20
 - "cross-key" 7, 19
 - dispersive pad 8, 19
 - "Er Or" displayed 18, 28
 - output 8, 20
- Extender electrode 35
- External filter 8, 17, 36
 - installing 23
- External lesion electrodes 35
- F**.....
- Faradic effect. See Electric shock
- Filters
 - external filter pack 23
 - installing 23
 - internal 23
 - specifications 33
 - warnings 4
- FINESSE filter pack 23
- FINESSE internal filter 23
 - changing 26
 - specifications 34
- Flow rate
 - vacuum 17, 33
- Footswitch 36
 - receptacle 21
- Frequencies
 - electrosurgical 16
- Frequency 38
- Front panel connectors 20
- Front panel controls 17
- Fulguration 10, 15, 38, 39
- Fuses
 - location 21
 - specifications 30
- G**.....
- Grounding pad. See Dispersive pad
- Guidelines
 - procedure 8
- H**.....
- Handswitch pen 36
- Handswitch receptacle 20
- Hemostatic techniques 10
- HEPA filter 16
- I**.....
- Indications for use 1
- Indicators
 - Error 19
- Inspection 1
- Inspection and preventive maintenance 26
- Interference
 - with other equipment 3
- Internal filter 8, 17, 34, 36
 - installing 23
- Interrupted sinusoid 39
- Isolated leads 38
- J**.....
- Jewelry
 - warnings 2
- L**.....
- Leads
 - warnings 2
- Leakage current 38
- LETZ
 - indications and contraindications 1, 10
 - post-procedure guidelines 13
 - procedure and technique 11
- LETZ electrodes 35
- Load resistance 39
- Loop electrode
 - selection 12
- loop electrodes 35
- M**.....
- Main power switch 20
- Maintenance
 - annual 26
 - cleaning 26
 - daily 26
- Maximum Output Voltage 3, 30, 39
- Metal objects
 - patient contact with 2
- Modulation 39
- Monitoring leads
 - warnings 2
- Monopolar electrosurgery 15, 39
- N**.....
- Needle electrodes 36
- O**.....
- Obtaining service 49
- Odor
 - adsorption 37
- Ohm 39
- Operation
 - normal operation of Finesse system 25
- Operational environment 34
- OptiMicro Needle electrodes 36
- Output control knob
 - coag 19
 - cut/blend 18
- Output error 8, 20
- Output settings 9
 - inadequate 3, 27
 - LETZ 10
- Output specifications 30
- Overcurrent protection 30
- P**.....
- Pacemakers
 - warnings 4
- Patient
 - contact with metal objects 2
- Power cord 21
- Power deficiencies
 - troubleshooting 27
 - warnings 3, 18, 19
- Preventive maintenance 26

R.....

- Radio frequency 39
- Rated Accessory Voltage 3, 32, 39
- Rear panel controls and connectors 21
- Receptacle
 - footswitch 21
 - handswitch 20
 - power cord 21
- Regulatory classifications 29
- REM. See CQM
- Repair
 - returning for 49
- Replacement parts and accessories 35

S.....

- Service 49
- Settings
 - output 9
- Setup 23
- Shutdown
 - post-procedure 25
- Sinusoid waveform 39

- Smoke evacuation system
 - description 8
 - principles 16
 - specifications 33
 - troubleshooting 27
 - Vacuum level switch 17
- Smoke filter connection 17
- Spark 39
- Specifications
 - audio tones 33
 - CQM circuit 32
 - environmental 34
 - fuses 30
 - output 30
 - physical dimensions 29
 - supply voltages and capacity requirements 29
- Spray coagulation 39
- Storage and transport environment 34
- Supplies
 - LETZ 11
- Supply voltages and capacity requirements 29

T.....

- Troubleshooting 27
 - error indications 19

U.....

- ULPA filter 16
- UtahBall electrodes 35
- UtahLoop electrodes 35

V.....

- Vacuum Level Switch 17
- Volt 39
- Volume control 21

W.....

- Warranty 49
- Watt 40
- Waveform
 - sinusoid 39
- Waveforms
 - electrosurgical 16

United States

Utah Medical Products, Inc.
7043 South 300 West
Midvale, Utah 84047
800.533.4984 • +1 801 566 1200

Ireland

Utah Medical Products Ltd.
Athlone Business & Technology Park
Dublin Road
Athlone, Co. Westmeath, Ireland
+353 (090) 64 73932

United Kingdom

Femcare-Nikomed Ltd.
Stuart Court, Spursholt Place, Salisbury Road
Romsey, Hampshire, SO51 6DJ
+44 (01794) 525 100

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Copyright © 2013 by Utah Medical Products, Inc.
FINESSE® , UtahLoop® UtahBall® , Safe-T-Gauge® , DXTender® , LETZ® , and Epitome® are registered trademarks, and FINESSE+™ , OptiMicro™ , and ZapGuard™ are trademarks, of Utah Medical Products, Inc., Midvale, Utah.