## **Test Procedures**

This checklist is intended to provide information for the preventive maintenance (PM) of the FINESSE+ Electrosurgical Generator and Smoke Evacuation Systems. It is intended for use by technical personnel who have experience with electrosurgical generator operation and maintenance, and who possess the test equipment and tools necessary to obtain the requested data. This document may not reference all tests necessary to demonstrate compliance to all regulations, nor should all tests listed be considered necessary for the safe use of the FINESSE+. Technical personnel should consult local regulations and industry recommendations to develop their own protocol for inspection and PM of the FINESSE+.

## **Required Equipment and Tools**

- Electrosurgical analyzer and leads
- Electrical safety analyzer, with patient lead connections
- Two-button electrosurgical switchpen and/or two-pedal footswitch
- SSE-500 FINESSE+ Internal Smoke Evacuation Filter

Mechani	ical Inspection and Service					
☐ Pass	Verify front panel controls, switches, and dispersive pad receptacle/pins are not damaged.					
☐ Pass	Verify switchpen receptacles firmly hold switchpen.					
	Knobs and color coded knob caps are securely attached.					
☐ Pass	Verify power cord and its rear panel receptacle are not damaged.					
☐ Pass	Verify footswitch cord (if present) and its rear panel receptacle are not damaged.					
$\square$ Done	Replace FINESSE+ Internal Smoke Evacuation Filter					
Performa	ance Inspection					
	Green LED in main power switch lights when switch is toggled to the "   " position.					
	One and only one cut mode indicator ("C", "1", "2", or "3") is illuminated					
	For all four cut/blend modes, the cut mode waveform selector knob adjusts so that light					
	indicator for each mode illuminates when knob's index pointer corresponds to the					
	correct indicator light.					
Verify po	wer displays can be adjusted through the full output range					
	Cut: 05-99					
	button on switchpen (repeat with cut pedal on footswitch, if available), and verify:					
□ Pass	Yellow cut activity indicator 😯 lights and cut mode tone is audible					
	Smoke evacuation motor activates, and remains running 5 seconds after deactivation					
	ag button on switchpen (repeat with coag pedal on footswitch, if available), and verify:					
☐ Pass	Blue coag activity indicator ights and coag mode tone is audible					
	Smoke evacuation motor activates, and remains running 5 seconds after deactivation					

	Pass Verify smoke evacuation flow rate switch changes motor speed.  Done Adjust volume control so that cut and coag audible tones can be clearly heard over the smoke evacuation motor set at the 'high' flow rate position.						er the	
<b>Dutput Calibration</b> Set cut mode output to "99", mode to "C", and electrosurgical analyzer load to $500\Omega$ . Line voltage MUST be $114.0-116.0$ VAC for FIN- $110$ and $228.0-232.0$ VAC for FIN- $220$ . Activate FINESSE+ system and adjust R334 until output is $97W$ . NOTE: Output may vary several watts due to feedback activity.								
□ Pass		(	(95W-99W)					
made thi more tha	rough a 500Ω	load. Tole for each re	rances are leading, and	isted in par allow at lea	rentheses. [	o not a	ements should to etivate system f een readings.	
□ Pass	Cut	30:	(25-35W)	60:	(51-69W)	90:	(77-103W)	
□ Pass	Blend 1	30:	(25-35W)	60:	(51-69W)	90:	(77-103W)	
□ Pass	Blend 2	30:	(25-35W)	60:	(51-69W)	90:	(77-103W)	
□ Pass	Blend 3	30:	(25-35W)	60:	(51-69W)	90:	(77-103W)	
□ Pass	Coag	30:	(25-35W)	60:	(51-69W)	75:	(64-86W)	
Set cut o connecti dispersiv	ng a non-indo	" and coag outive 200Ω open. Activ	2 load betw	een the act	ive electrod	e and m	urrent path by ains ground. Le eakage current:	
	e HF Leakage Cut				Pass Coag		(<150mA)	
Connect non-inductive $200\Omega$ load resistor between dispersive electrode and mains ground. Leave active connection open. Activate FINESSE+ system to check dispersive HF leakage current: Dispersive HF Leakage								
	Cut	_			Pass Coag		(<150mA)	

Touch C	urrents					
listed in	electrical safety analyze parentheses: nd lead intact	r, measure the follo	owing con	ditions. A	.cceptable cond	ditions are
□ Pass Grou	Normal Polarity nd lead open					
☐ Pass	Normal Polarity	_(<100μA)	☐ Pass	Reverse F	olarity	_ (<100µA)
appropria	e electrical safety analyz ate connections for the f g the FINESSE+ system.	ollowing leakage cu	urrents. N	/leasure th	ie leakage curr	ents without
☐ Pass	Patient Auxiliary Curren	t	(<100µA)	)		
☐ Pass	Patient Leakage Current	: (PLC)	(<100µA)	)		
□ Pass	PLC Caused by Ext Volta (aka Mains on Applied F	•	(<5000µ/	A)		
	rcuit Tests					
	ctivate the FINESSE+ sys	-	sts, as mo	st ES analy	zers are not de	esigned to
	d activation during CQM g the instructions for the	•	rm the fol	lowing me	easurements IN	I THE ORDER
	vithout disconnecting th	•		_		
resistor s	ubstitution box if analyz	er does not have C	QM test c		•	
	INESSE+ system and tur	•	witch.			
	Verify $\uparrow \Box$ error is active					
☐ Pass	CQM dispersive pad illu	minates the CQM p	ad icon $\widehat{}$	<b>a</b> bov	e the pad rece	ptacle.
☐ Pass	Reduce resistance until	nt error clears.		130	)Ω(11	7Ω-143Ω)
☐ Pass	Reduce to $50\Omega$ , then ra	ise until 🔽 error	activates.	689	Ω (61	$\Omega$ -75 $\Omega$ )
□ Pass	Reduce resistance until	error activates	5.	100	$\Omega$ (5 $\Omega$	2-15Ω)
	ror Checks					
☐ Pass	Pressing both cut and c	-	ons or foo	otswitch p	edals simultan	eously
	produces error and				. 4	
⊔ Pass	Activation of coag mod	e into open circuit	conditions	s does not	produce 📆	error

Notes: Performed by: FINESSE+ Serial No: Signature: Date:

Test Equipment	Make/Model	Serial/Reference	Calibration Expiry
ES Analyzer			
Safety Analyzer			